

Factors Affecting Orthoses Selection:



Jerome Recommended Cervical Orthotics

Location	Clinical Condition	Miami J	JTO	V ₁ Halo	Surgery
All	Cervical Strain	✓			
High (C1-C2)	C1 Ring (Jefferson Fx)				
	Stable	✓			
	Unstable			✓	
	Odontoid Fx Type I		✓		
	Type II & III			✓	
	C2 Neural Arch (Hangman's Fx)			✓	
Mid (C3-C5)	Flexion Injuries	✓			✓
	Extension Injuries	✓		✓	
Low (C5-T2)	Flexion Injuries		✓		✓
	Extension Injuries		✓		✓