


Universal Halo Adapter

Meet the challenge of intra-operative stabilization with Jerome's **Universal Halo Adapter**.

Part No. 581

The UHA bridges the gap between the Mayfield™ Neurosurgical Headrest Arm and any halo ring.

Safely and securely position the patient for anterior and posterior surgical approaches and reduce the risk of unwanted patient movement.



*Anterior Approach
with ReSolve® Halo Ring*

Posterior Approach



- Versatility to adapt to both Anterior and Posterior Approaches.
- Compatible with all halo rings.
- Halo Attachment Pin inverts to fit standard Titanium Pins and ReSolve® pins.
- Streamlined design does not encroach the surgical field and permits full access to patient for anesthesia.

www.jeromemedical.com

305 Harper Dr., Moorestown NJ USA 08057-3239 • 800-257-8440 • Fax: 856-778-5333

 **JEROME**
CERVICAL SPINE SYSTEM®
Helping You Meet The Challenge